

### Reuse and reprocess of surgical saw blades

Disposable product marked ② may not be reused or reprocessed, these products are designed for single use only.

Reusing and/or reprocessing negatively affects the performance of these devices. All material is held to exacting tolerances, which is not guaranteed if reused. The teeth, body, and hub are designed to specification, and dimension will be irreversibly changed during use.

Reusing and/or reprocessing may cause damage to the surface of the product, leading to corrosion of the device. To protect patients, O.R staff, and surgeons from infection of contaminated devices, products are single use only.

Product is sterile only if original packaging is unopened and undamaged. If the packaging is opened or damaged, the device is no longer sterile and should be returned to the supplier.