

Reuse and reprocess of eLIBRA[®] Soft Tissue Force Sensor Devices

Disposable product marked ② may not be reused or reprocessed, these products are designed for single use only. eLIBRA[®] Soft Tissue Force Sensors are factory calibrated devices. Any reuse and/or reprocessing will cause damage to the device rendering it out of calibration and/or specification.

To protect patients, OR staff, and surgeons from infection of contaminated devices, products are designed for single use only. Product is sterile only if original packaging is unopened and undamaged. If the packaging is opened or damaged, the device is no longer sterile and should be returned to the supplier.